

Bringing flexibility and affordability to CPAP therapy with a truly compact and portable CPAP system from evo. This CPAP has a pressure from 4-20 cmH₂O, patient-adjustable ramp, automatic altitude adjustment, leak compensation, international power supply, compliance meter, simple set-up and the optional DataRest™ Compliance System.

R_x ONLY



Compact & Portable

everest[®] 2 Includes

- CPAP
- 22 mm hose, 6'
- Power supply and power cord
- 2 pack air filters
- Custom system travel bag
- Interface Pouch
- Patient Manual

everest[®]2 Accessories

Re-order

802CH	CPAP & Heated Humidifier
802H	Heated Humidifier
802B	Battery
802MPA	Mobile Power Adapter
802S	DataRest Software Kit
802PK	DataRest Patient Kit
802DC	DataRest Data Card

everest[®]2 Model Numbers

EV2	CPAP
HH1	Heated Humidifier
E2H	CPAP & Heated Humidifier System

Device Specifications

Pressure Range:	4-20 cmH ₂ O
Input Power:	100-240 VAC, 50-60Hz, 12 VDC
Weight:	2.1 lbs (1 kg)
Size:	5"x6"x5"
Warranty:	2-Year CPAP 2-Year Heated Humidifier 1-Year Battery 3-Year Mobile Power Adapter



evo Medical Solutions

26378 289th Place
Adel, IA 50003
515.993.5001
515.993.4172 f
800.759.3038 (toll free)