



# APAP

## APAP & Humidifier System

804

Our APAP machine joins the sleep product line — expanding choices for our sleep apnea patients. The 804 APAP monitors the user's breathing performance with pressure sensors and works when the user needs it. The APAP is considered one of the most advanced forms of treatment for sleep apnea available today.

R<sub>x</sub> ONLY



APAP & Humidifier

### APAP Includes

804 APAP & Heated Humidifier, detachable power cord, flexible tubing, 2 pack filters, mini USB port, mini SD card port, carry bag and instruction manuals

### Device Specifications

Pressure:	4-20 cmH2O, ≤30 DBA @10 cmH2O, 1 meter
Ramp Time:	0-45 min
Ramp Starting Pressure:	3-19 cmH2O
APAP Size:	5.8"x5"x3.9"
Heater Size:	10"x6.8"x1.6"
Water Chamber:	5.9"x4.5"x4.2"
APAP Weight:	1.76 lbs.
Heater Weight:	.84 lbs.
Water Chamber:	.53 lbs.
Power Supply:	Universal power supply 100V-240V, 50/60Hz
Humidifier Power Supply:	Universal type 120V-240V, 50/60Hz
Humidifier Water Capacity:	450 ml
Warranty:	2 years APAP 1 year Humidifier



### evo Medical Solutions

26378 289th Place  
Adel, IA 50003  
515.993.5001  
515.993.4172 f  
800.759.3038 (toll free)  
www.evomedical.com