

Disposable Canister

800cc Disposable Suction Canister

610

This superior 800cc Disposable Suction Canister is available in either 12 or 48 count cases. This superior 800cc Disposable Suction Canister has a unique float shut off that helps prevent fluid from entering the aspirator.

Collection Canisters

610-12 800cc Disposable Suction Canister (12/case)

610-48 800cc Disposable Suction Canister (48/case)

Collection Canister With Accessories

22330 Canister Kit (each)



Canister Kit



evo Medical Solutions

26378 289th Place

Adel, IA 50003

515.993.5001

515.993.4172 f

800.759.3038 (toll free)