

MediAspirator

Portable Suction Unit

606

The MediAspirator 606 is the ideal suction unit designed for inpatient and outpatient care requirements. One collection jar and one connection tubing are included.

Features:

- high-performance vacuum specifications for most home care suction applications
- durable piston type compressor
- easy-to-read gauges

R_x ONLY



800cc Disposable Canister



Canister Kit

MediAspirator Includes

800cc Disposable Suction Canister
6' Connecting Tubing
Hydrophobic filter
Instruction manual

MediAspirator Accessories

T-12 Collection Jar Tubing package of 10
612 Complete Air Filter Assembly
613 6' Aspirator Tubing
610E Right Angle Connector
611 Hydrophobic filter
610-12B 800cc Disposable Suction Canister (12/case)
610-48B 800cc Disposable Suction Canister (48/case)
22330 Canister Kit (each)

Device Specifications

Dimensions (LxHxW in): 13.4" L x 8.3" H x 6.1" W
(34cm L x 21cm H x 15.5cm W)

Weight: 9.5lbs. (4.3 kg)

Secretion container volumes: filling limit: 27oz. (800 ml)

Max. vacuum: 23.6" Hg~26.4" Hg (600~700 mmHg)

Noise level: $\leq 50\text{Db}$

Power supply: 120V/60 Hz, 230v/50 Hz

Flow: $\leq 16\text{LPM}$ Maximum

Warranty: 1-Year Device



evo Medical Solutions

26378 289th Place

Adel, IA 50003

515.993.5001

515.993.4172 f

800.759.3038 (toll free)