

# RemRest

Heated Humidifier

910

The RemRest Heated Humidifier integrates directly with the evo RemRest 900 series CPAP devices. It has a removable water reservoir that is reusable and easy to clean. The humidifier moistens and warms the air delivered by the CPAP devices.

R<sub>x</sub> ONLY



## RemRest Heated Humidifier Includes

Water Reservoir  
Extra Frame Support  
Manual

## Accessories

911 Water Reservoir

## Device Specifications

Size:	8.4" W x 13.1" L x 4.7" H
Weight:	2.5 lbs.
Water Capacity	1 <sup>2</sup> / <sub>3</sub> cups (400 ml) at maximum water level
Warranty:	1 Year on device
Heater Settings: (Plate Temp.)	1-5 (85°F to 149°F) (29.4°C to 65°C)
Humidity Output:	30 g H <sub>2</sub> O/hr @ 25 lpm w/ambient air @ 72°F 20% RH
Pressure Drop with Humidifier:	0.1 cm H <sub>2</sub> O @ 40 lpm flow