

Nebulizer

Reusable Jet Nebulizer Set

22226

Product Features:

- 6 cc medication cup with graduations
- Compatible with standard masks
- Billing code: A7005
- Dishwasher safe
- Latex Free

R_x ONLY

Includes

Reusable Nebulizer
Tubing
Mouthpiece
T-Connector
Reservoir Tube
Instructions for use

Device Performance

High nebulation rate: (3cc @ 7 min.)
>80% of particles below 5 microns

Warranty: 6 months



evo Medical Solutions

26378 289th Place

Adel, IA 50003

515.993.5001

515.993.4172 f

800.759.3038 (toll free)