

TranSport

Portable Compressor Nebulizer

5055

The TranSport is a compact, piston compressor providing high-performance while weighing only 370g. Versatility. Dependability. Simplicity. Convenience.

R_x ONLY



TranSport with Reusable Cup



Rechargeable Battery



Carry Bag



Auto Adapter

TranSport Includes

5055 TranSport Compressor/Nebulizer
Reusable Medication Cup
Carry Bag
Power Cord
Operating Manual

TranSport Accessories

5055BAT	Rechargeable Battery	5055MA	Mask Adapter
5055CORD	Auto Adapter	5055ADP	Adapter
5055BAG	Carry Bag	5055MP	Mouthpiece
5055AM	Adult Mask		
5055PM	Pediatric Mask		
5055A	Tubing Adapter		
5055F	Replacement Filters (5pk)		
5055MC	Reusable Medication Cup		

Device Specifications

Maximum Pressure:	≥30 PSI
Operating Flow Rate:	≥3.5 lpm @ 10 PSI operating pressure
Warranty:	3 Year device 180 day Battery
Max Dosage:	6cc
Weight (compressor only)	.826lbs/1.5lbs (370g/675g) w/battery

TranSport Battery Specifications

Battery Type:	NIMH 1500 MAH
Supply Voltage:	15VDC Use only original TranSport Adapter
Complete Charge:	At least 4 hours
Battery Life:	Average of two 20 minute treatments at full charge



evo Medical Solutions

26378 289th Place
Adel, IA 50003
515.993.5001
515.993.4172 f
800.759.3038 (toll free)