

VacuMax™

The Ultimate Portable Aspirator

605

The portable VacuMax weighs only 7 lbs and comes complete with its own carry bag. This unit operates on AC, DC or internal battery, making it ideal for patient mobility or as back-up to a larger AC only unit.

R_x ONLY



800cc Disposable
Canister



Canister Kit

VacuMax 605 Includes

- Carry Bag
- 800cc Disposable Suction Canister
- 6' Aspirator Tubing
- 7" Collection Jar Tubing
- 5 Unit Filters
- Hydrophobic filter
- Instruction manual
- DC power adapter
- AC Power Cord
- Internal battery

VacuMax Accessories

- 605A VacuMax Carry Bag
- 20289 Auto Adapter
- 612 Complete Air Filter Assembly
- 613 6' Aspirator Tubing
- T-12 7" Collection Jar Tubing (10 pk)
- 1005 External Filter (20 pk)
- 610E Right Angle Connector (20 pk)
- 611 Hydrophobic filter
- 610-12 800cc Disposable Suction Canister (12/case)
- 610-48 800cc Disposable Suction Canister (48/case)
- 22330 Canister Kit (each)

Device Specifications

Suction:	20 inHg (50.8 cm Hg) minimum
Flow:	12 LPM minimum with Hydrophobic Filter
Weight:	7 lbs 3 oz (3.2kg)
Warranty:	1-Year device 90 day battery



evo Medical Solutions

26378 289th Place
Adel, IA 50003
515.993.5001
515.993.4172 f
800.759.3038 (toll free)