

Ventlab Nasal CPAP Cannula System consists of headgear, padding, circuit, tubing clips and a single nasal cannula. Two pieces of flexible tubing join the cannula to a "Y" shaped coupling for connection to a positive pressure device. The system includes ports for pressure monitoring, gas sampling and O2 delivery.

Features

- Latex Free
- Soft slimline circuit
- Easy to clean
- Wider field of vision
- Comfort fit nasal inserts
- Multiple adjustments for greater stability and comfort
- Comfortable featherlite padded headgear

Ventlab Nasal CPAP Cannula System contains

- 4 sizes of nasal inserts
- 2 sizes of adjustable headgear
- Patient circuit with swivel connector
- Tubing clips
- Headgear padding
- Convenient nylon carrying case

Specifications

Pressure vs. Flow

1.3 cm H₂O @ 40 lpm

2.6 cm H₂O @ 60 lpm

5.4 cm H₂O @ 80 lpm

6.2 cm H₂O @ 100 lpm

Component Re-ordering Information

CS9000	Nasal CPAP System
CS9001	Large Headgear w/clips and padding
CS9002	Small Headgear w/clips and padding
CS9003	Small Cannula
CS9004	Medium Cannula
CS9005	Large Cannula
CS9006	X-Large Cannula
CS9007	Patient Circuit
CS9008	Tubing Clips
CS9009	Headgear Padding
CS9011	Soft Headgear



evo Medical Solutions

26378 289th Place

Adel, IA 50003

515.993.5001

515.993.4172 f

800.759.3038 (toll free)