

VixOne™

High Performance Disposable Nebulizer

A study comparing the performance of various disposable nebulizers concluded that the VixOne provided the largest respirable mass of aerosol for the patient.

R_x ONLY

Ordering Information

21304	VixOne Nebulizer, Case of 50
21306	VixOne Nebulizer, w/Adult Mask, Case of 50
21308	VixOne Nebulizer, w/Pediatric Mask, Case of 50
1010	VixOne Nebulizer, w/Dragon Mask, Sold individually
1007	Pediatric Mask
1008	Adult Mask

Features

10µm particle size production (MMAD = 2.7 micron) allowing high lung desposition which results in maximum therapeutic effect.*

- Anti-spill "tee" design
- Leak proof and secure seal with simple 1/4 turn assembly
- Low residual volume (Ø.7cc)
- High medication bowl capacity (13cc)
- Nebulizes in upright position or at angles up to 45°

* Reference (Hest Vol. 110, p.498-505 August 1996)




VixOne w/Dragon Mask

Nebulizer Output vs. Mass Available to Patient

■ Respirable Mass Output (mg) ■ Respirable Mass Available (mg)

Dean Hess et al CHEST/110/August 1996



evo Medical Solutions

26378 289th Place
Adel, IA 50003
515.993.5001
515.993.4172 f
800.759.3038 (toll free)